

Sample Name :Glucagon (1-29), Human
 Cat No: RP10772
 Lot No : E4934720L
 Time Processed :11:59:30
 Month-Day-Year Processed :12/24/2020

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	15
25.00	Pumps	B.Conc	75
25.01	Pumps	B.Conc	95
31.00	Pumps	B.Conc	95
31.01	Pumps	B.Conc	15
40.00	Pumps	B.Conc	15
40.00	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
 Equipment: ZJ19010140

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.900	255823	28669	0.921
2	9.738	151257	25673	0.544
3	10.022	26798396	1823289	96.463
4	10.515	69319	8538	0.250
5	10.645	47621	5911	0.171
6	11.070	93210	10173	0.336
7	12.154	365404	45997	1.315
Total		27781031	1948250	100.000