


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 05/10/18  
Time Processed : 10:55:20  
Injection Volume : 0.2  
Sample Name : Rabies Virus Glycoprotein  
Sample ID : U8495DE070-1  
Theoretical MW : 3266.70  
Observed MW : 3267.00

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H2O/50% MeOH