


Mass Spectrum



Sample Information

Month-Day Processed : 08/20/20
Time Processed : 08:05:03 PM
Injection Volume : 0.3
Sample Name : Jagged-1 (188-204)
Cat No. : RP20331
Lot No. : E6903820H
Theoretical MW : 2107.40
Observed MW : 2106.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50% H_2O /50% MeOH