


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/15/19
Time Processed : 06:00:10 PM
Injection Volume : 0.3
Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. : RP11110
Lot No. : E6805719J
Theoretical MW : 2197.60
Observed MW : 2197.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH