


# Mass Spectrum



Sample Information  
Month-Day Processed : 07/08/21  
Time Processed : 21:40:59  
Injection Volume : 0.4  
Sample Name : LCMV gp33 - 41  
Cat No. : RP20257  
Lot No. : E6927921G  
Theoretical MW : 1044.25  
Observed MW : 1044.0

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH