


Mass Spectrum



Sample Information
Month-Day Processed : 01/16/21
Time Processed : 8:19:44 AM
Injection Volume : 0.3
Sample Name : Exendin-4
Cat No. : RP10874
Lot No. : E7242421A
Theoretical MW : 4186.66
Observed MW : 4185.9

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH