


Mass Spectrum



Sample Information
User : Gary
Month-Day Processed : 08/29/19
Time Processed : 19:08:18
Injection Volume : 0.3
Sample Name : Calcitonin Gene Related Peptide (CGRP) (8-37), rat
Cat No. : RP11090
Lot No. : E4087419H
Theoretical MW : 3127.60
Observed MW : 3127.2

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200