


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/14/18
Time Processed : 16:32:23
Injection Volume : 0.2
Sample Name : TAT peptide
Cat No. : RP20256
Lot No. : PE2531803
Theoretical MW : 1621.91
Observed MW : 1622.10

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH