


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 09/27/19  
Time Processed : 1:13:14 PM  
Injection Volume : 0.8  
Sample Name :  $\beta$ -Amyloid(1-42),human  
Cat No. : RP20527  
Lot No. : E1274319I  
Theoretical MW : 4514.10  
Observed MW : 4514.0

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250 $\pm$  $\epsilon$   
Block Temp : 200 $\pm$  $\epsilon$

Equipment: ZJ17010408  
Interface Bias: +4.5kV  
Drying Gas Flow: 5 L/min  
T.Flow: 0.2ml/min  
B.conc: 50%H<sub>2</sub>O/50%MeOH