


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 06/25/20
Time Processed : 12:01:21 PM
Injection Volume : 0.4
Sample Name : Myelin Oligodendrocyte Glycoprotein Peptide (35-55), mouse, rat
Cat No. : RP10245
Lot No. : E5578720F
Theoretical MW : 2581.98
Observed MW : 2581.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH