


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 05/19/17
Time Processed : 15:01:50
Injection Volume : 0.2
Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. : RP11110
Lot No. : PE6691705
Theoretical MW : 2197.60
Observed MW : 2197.50

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH