


Mass Spectrum



Sample Information

Month-Day Processed : 05/23/21
Time Processed : 14:31:29
Injection Volume : 0.2
Sample Name : MDYKDHGQDYKDHIDYKDDDDK peptide
Cat. No. : RP21087
Lot. No. : E6789321E
Theoretical MW : 2861.87
Observed MW : 2861.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH