


# Mass Spectrum



## Sample Information

Acquired by :	Gary	Equipment	: GK11010007
Month-Day Processed :	03/30/19	Interface Bias	: +4.5 kV
Time Processed :	01:47:50 PM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.2	T.Flow	: 0.2 ml/min
Sample Name :	OVA Peptide (323-339)	B.conc	: 50% H2O/50% MeOH
Cat No :	RP10610		
Lot No :	PE5791903		
Theoretical MW :	1773.91		
Observed MW :	1773.9		