


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 05/18/17
Time Processed : 7:47:05
Injection Volume : 0.2
Sample Name : c-Myc Peptide
Cat No. : RP11731
Lot No. : PE4171705
Theoretical MW : 1203.30
Observed MW : 1203.50

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH