


Mass Spectrum



Sample Information
Month-Day Processed : 07/02/21
Time Processed : 11:54:28
Injection Volume : 0.2
Sample Name : Exendin-4
Cat No. : RP10874
Lot No. : E8894121F
Theoretical MW : 4186.66
Observed MW : 4186.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH