

Certificate of Analysis

Product Information

Product Name: PIVKA II Antibody (2D7), mAb, Mouse
Cat. No: V06702
Quantity: Specified on product label
Presentation: 50 mM Na-citrate, 150 mM NaCl, pH7.0, containing 0.03% Proclin 300.
Storage: For up to three months at 2-8°C or for up to three years at -20°C or below. Avoid repeated freezing and thawing cycles.
Reactivity: Specific to human PIVKA II
Purification: Protein A/G affinity column
Lot No.: T1902001

QC Results

Characteristics	Methods	Limits	Values
Concentration:	A280	≥ 2 mg/mL	5.8 mg/mL
Purity:	SDS-PAGE,SEC-HPLC	≥ 95%	98.3%
Aggregation:	SEC-HPLC	≤ 5%	1.7%
Isoelectric Point:	cIEF	match	N/A
Expiry Date:			02/22/2022

Date : 02/22/2019

Approved by: 

Antibody Production Director

For research use only